


Amplatzer™ Septal Occluder

OPTIMAL OCCLUSION, **PROVEN RESULTS**



STRUCTURAL HEART | STRUCTURAL INTERVENTIONS



WHAT HAS MADE THE AMPLATZER™ SEPTAL OCCLUDER (ASO) **THE WORLD'S LEADING AND MOST TRUSTED** DEVICE FOR ATRIAL SEPTAL DEFECT (ASD) CLOSURE FOR OVER 25 YEARS?

1 PROVEN SAFE AND EFFECTIVE CLOSURE

Since ASO's first formal approval in 1998, over 500,000 patients¹ across the globe have benefitted from this innovative technology.

ASO has long been considered the standard of care, with a clinically validated 98.5% closure rate² and extremely low major and minor complications.

2 CLINICAL EVIDENCE AND EXPERTISE YOU CAN TRUST

No other ASD closure device comes close to having the amount of clinical data gleaned over the past quarter-century as ASO. In addition to multiple clinical trials, there have also been over 750 peer-reviewed articles about its safety and effectiveness. By contrast, some competitive devices have no clinical trial data of their own, and in some cases have leveraged ASO's data for their regulatory approvals.


With nearly 1,400 combined years of ASD closure experience, our product support reps are trusted partner of interventional cardiologists who provide unmatched case support and education.

3 SMART PRODUCT DESIGN AND EASE-OF-USE

The braiding pattern of the ASO has an optimal design for defect closure as well as allowing access for future left side heart interventions.

The ASO device is compatible with the Amplatzer™ Trevisio™ Intravascular Delivery System, providing a simple procedure and confident septal alignment prior to device release.

ASO offers the widest range of sizes on the market – from 4 mm up to 38 mm³ – so you can find the right device for every patient.




1 PROVEN SAFE AND EFFECTIVE CLOSURE

Backed by PIVOTAL study among others, plus 750+ peer reviewed journals

Since its introduction over 25 years ago, ASO has consistently provided safe and effective closure, with minimal adverse events^{2, 4}

NOTE: Comparing to competitive products is often challenging due to their lack of independent clinical studies.

EASE-OF-USE AND HIGH PROCEDURAL SUCCESS AS SHOWN IN THE US PIVOTAL TRIAL 1998-2001²



ASO LEADS IN SAFETY WITH A SUPERB CLINICALLY-VALIDATED SAFETY PROFILE AND MINIMAL ADVERSE EVENTS

EROSION RATE³

0.043%-0.3%

No erosions reported for sizes smaller than 13 MM. Erosion can be mitigated by sizing technique.

DEVICE EMBOLIZATION² with percutaneous removal

0.2%

NOTE: A competitive ASD product shows a 20x higher embolization rate.⁵

WIRE FRAME FRACTURES²

0.0%

NOTE: A competitive ASD product study at 36 months shows wire frame fracture rate of 57% with rates ranging from 38% to 89.5% depending on the device size.⁶

NEW ARRHYTHMIAS² requiring major treatment

0.5%

NOTE: A competitive ASD product showed a real-world rate of new arrhythmias requiring hospitalization or intervention of 4.8%.⁷

2 CLINICAL EVIDENCE AND EXPERTISE YOU CAN TRUST

No ASD device in the world has been more studied than ASO

Unmatched clinical evidence

The Amplatzer™ Septal Occluder is the most studied transcatheter atrial septal defect (ASD) closure device available today, with over 25 years of demonstrated clinical experience.

THESE 3 CLINICAL STUDIES ENCOMPASS OVER 1,900 ASO PATIENTS

- U.S. PIVOTAL Trial: 442 patients
- Amplatzer Septal Occluder Post-Approval Study: 1000 patients
- Magic Atrial Septal Defect Study: 478 patients



AMPLATZER™ SEPTAL OCCLUDER DEVICES SOLD¹



PATIENTS STUDIED^{3, 4, 8}



PEER-REVIEWED ARTICLES

UNMATCHED PRODUCT AND INDUSTRY EXPERTISE



WORLD-CLASS EDUCATION AND SUPPORT

Providing expert instruction to ensure your team has convenient access to the latest training



LARGEST PORTFOLIO

Offering the broadest product portfolio focused on helping people live fuller, healthier lives



EMPATHIC PATIENT SUPPORT

Providing world-class patient materials and resources to ensure positive experiences



CLINICAL TRIALS

Ongoing world-class research to balance risk/benefit ratio, improve safety and advance standard-of-care

3 SMART PRODUCT DESIGN AND EASE-OF-USE

with key product features that support excellent outcomes

TREATMENT RANGE 3 MM TO 38 MM


No other ASD closure device on the market today can match ASO's wide range of device sizes.

AVAILABLE IN 26 DIFFERENT SIZES³

Ranging from 4 mm to 38 mm, designed to treat most patient anatomies and constructed to accommodate various septal morphologies.



UNIQUE PRODUCT FEATURES




DOUBLE DISCS are designed to securely appose both sides of the septal wall.

The **WIDE CONNECTING WAIST** is designed to center the device and fill the ASD.

SHAPE-MEMORY NITINOL MESH
Designed to securely appose both sides of the septal wall.

INTAGLIO™
A proprietary wire treatment for Amplatzer™ Nitinol wire. The only treatment proven to reduce nickel leaching by more than 95%¹.

POLYESTER MATERIAL
Promotes occlusion and tissue in-growth.



Enables left atrium access for future procedures including ablation procedures, mitral valve procedures, and left atrial appendage closure.⁹⁻¹¹

AMPLATZER™ TREVISIO™ INTRAVASCULAR DELIVERY SYSTEM

From delivery to deployment, physicians appreciate the Amplatzer Trevisio Delivery System's ideal combination of torque strength, sheath diameter, reliable performance and ease-of-use.

The Trevisio Delivery System is an ultra-flexible delivery system enabling interventional cardiologists to perform transcatheter ASD closures work with complete confidence.

Compatible with all Amplatzer Septal Occluders

- A unique one-piece cable design improves control, flexibility, and predictability
- Ideal for a wide range of heart morphologies

ULTRA-FLEXIBLE TIP
Improves assessment of device position before cable release and reduces bias on the device.

FLEXIBLE TRANSITION SECTION


Provides structure, pushability, and maintains sheath position during deployment of the device.

STIFF PROXIMAL SECTION

Provides reliable and steady pushability and control



WATCH THE TREVISIO DELIVERY SYSTEM IN ACTION



REFERENCES

1. Data on file at Abbott
2. Du, Hijazi, Kleinman, Silverman, Larntz. Comparison between transcatheter and surgical closure of secundum atrial septal defect in children and adults. *Journal of the American College of Cardiology*. 2002;39(11):1836-1844.
3. ARTEN600038247A ASO IFU
4. Turner DR, Owada CY, Sang Jr. CJ, Khan M, Lim S. Closure of Secundum Atrial Septal Defects With the AMPLATZER Septal Occluder: A Prospective, Multicenter, Post-Approval Study. *Circulation: Cardiovascular Interventions*. 2017;10(8):1-7.
5. Haas NA, Soetemann DB, Ates I, et al. Closure of Secundum Atrial Septal Defects by Using the Occlutech Occluder Devices in More Than 1300 Patients: The IRFACODE Project: A Retrospective Case Series. *Catheterization and cardiovascular interventions : official journal of the Society for Cardiac Angiography & Interventions*. 2016;88(4):571-581. doi:10.1002/ccd.26497
6. CSI America presentation by Dr. Qureshi, A. ASD Closure: Update on Cardioform Devices: CSI America, Nov 15, 2023. <https://www.csi-congress.org/conferences-courses/conferences/csi-america/archive>
7. GORE™ CARDIOFORM™ ASD Occluder Instructions for Use Feb 2023 IFU 30259376
8. Everett AD, Jennings J, Sibinga E, et al. Community use of the Amplatzer atrial septal defect occluder: results of the multicenter MAGIC atrial septal defect study. *Pediatr Cardiol*. 2009;30:240-247.
9. Santangeli P, Di Biase L, Burkhardt JD, et al. Transseptal access and atrial fibrillation ablation guided by intracardiac echocardiography in patients with atrial septal closure devices. *Heart rhythm*. 2011;8(11):1669-1675.
10. Villablanca PA, Lee J, Wang DD, et al. Transseptal Puncture Through an Amplatzer Atrial Septal Occluder for Edge-to-Edge Repair With MitraClip NTr System. *Cardiovascular revascularization medicine : including molecular interventions*. 2020;21(11S):63-64.
11. Gloekler S, Shakir S, Meier B. Transseptal Puncture Through Amplatzer Atrial Septal Occluder for Left Atrial Appendage Closure. *JACC Cardiovascular interventions*. 2017;10(21):2222-2223.

IMPORTANT SAFETY INFORMATION AMPLATZER™ SEPTAL OCCLUDER

R ONLY INDICATION FOR USE

The Amplatzer™ Septal Occluder is a percutaneous, transcatheter, atrial septal defect closure device intended for the occlusion of atrial septal defects (ASD) in secundum position or patients who have undergone a fenestrated Fontan procedure and who now require closure of the fenestration. Patients indicated for ASD closure have echocardiographic evidence of ostium secundum atrial septal defect and clinical evidence of right ventricular volume overload (such as, 1.5:1 degree of left-to-right shunt or RV enlargement).

CONTRAINDICATIONS

The Amplatzer™ Septal Occluder is contraindicated for the following: Any patient known to have extensive congenital cardiac anomaly which can only be adequately repaired by way of cardiac surgery; Any patient known to have sepsis within 1 month prior to implantation, or any systemic infection that cannot be successfully treated prior to device placement; Any patient known to have a bleeding disorder, untreated ulcer, or any other contraindications to aspirin therapy, unless another antiplatelet agent can be administered

for 6 months; Any patient known to have a demonstrated intracardiac thrombi on echocardiography (especially left atrial or left atrial appendage thrombi); Any patient whose size (such as, too small for transesophageal echocardiography probe, catheter size) or condition (active infection, etc.) would cause the patient to be a poor candidate for cardiac catheterization; Any patient where the margins of the defect are less than 5 mm to the coronary sinus, inferior vena cava rim, AV valves, or right upper lobe pulmonary vein.

POTENTIAL ADVERSE EVENTS

Potential adverse events may occur during or after a procedure placing this device may include, but are not limited to: Air embolus; Allergic dye reaction; Anesthesia reactions; Apnea; Arrhythmia; Cardiac tamponade; Death; Embolization; Fever Hypertension/hypotension; Infection including endocarditis; Need for surgery; Pericardial effusion; Perforation of vessel or myocardium; Pseudoaneurysm including blood loss requiring transfusion; Stroke; Tissue erosion; Thrombus formation on discs; Valvular regurgitation.

IMPORTANT SAFETY INFORMATION AMPLATZER™ TREVISIO™ INTRAVASCULAR DELIVERY SYSTEM

R ONLY INDICATIONS AND USAGE

The Amplatzer™ Trevisio™ Intravascular Delivery System is intended to provide a pathway through which devices are introduced within the chambers and coronary vasculature of the heart or in the peripheral vasculature.

CONTRAINDICATIONS

None known.

WARNINGS

- This device was sterilized with ethylene oxide and is for single use only. Do not reuse or resterilize this device. Attempts to resterilize this device can cause a malfunction, insufficient sterilization, or harm to the patient.
- Do not use this device if the sterile package is open or damaged.
- Use on or before the last day of the expiration month that is printed on the product packaging label.
- The sheath is designed to be used with the loader. Do not attach a syringe directly to the sheath because the sizing is incompatible and may result in ingress of air or excessive bleeding.
- Use the hemostasis valve to impede the backflow of blood during the implant procedure.
- Do not use a power injection syringe to inject contrast solution through the sheath.
- Remove the dilator and sheath from the patient slowly to prevent an ingress of air.

PRECAUTIONS

- Store in a dry place.
- This device should be used only by physicians who are trained in standard transcatheter techniques. The physician should determine which patients are candidates for procedures that use this device.
- The physician should exercise clinical judgment in situations that involve the use of anticoagulants or antiplatelet drugs before, during, and/or after the use of this delivery system.
- Use caution when advancing the dilator and sheath to avoid damaging tissue and vessels or interfering with previously implanted medical devices.
- Use standard transcatheter techniques when using Amplatzer™ products.

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur during or after a procedure using this delivery system may include, but are not limited to: Air embolism, Arrhythmia, Arteriovenous fistulae, Bleeding, Brachial plexus injury, Cardiac tamponade, Death, Dissection, Endocarditis, Hematoma, Infection, Myocardial infarction, Perforation, Peripheral embolism, Peripheral pulse loss, Stroke, Thrombosis, Tissue trauma/damage, Valve damage, Vascular occlusion, Vessel trauma/damage

CAUTION: Product(s) intended for use by or under the direction of a physician. Prior to use, reference to the Instructions for Use, inside the product carton (when available) or at <https://www.eifu.abbott/> for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

See Important Safety Information referenced within.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott.

Abbott

3200 Lakeside Dr., Santa Clara, CA. 95054 USA, Tel: 1.800.227.9902

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

www.structuralheart.abbott

© 2024 Abbott. All Rights Reserved. MAT-2400576 v3.0 | Item approved for U.S. use only.

